

KneeStim

Part No: K-KNEE

Conductive Garment for use with TensCare TENS & EMS units.

For single patient use only.

Universal size supplied with one pack of four reference electrode pads (50 x 100 mm electrodes).

TensCare Ltd

PainAway House,
9 Blenheim Road,
Longmead Business Park,
Epsom, Surrey KT19 9BE UK.
Tel: +44 (0)13 72 72 34 34
Web: www.tenscare.co.uk



Instructions for Use

1. Ensure that the stimulator is turned **OFF**.
2. Insert one pin connector (from one lead pair) into the connector attached to the KneeStim and the other (from the same lead pair) into the connector on the reference electrode pad.
3. The KneeStim should be worn slightly damp to improve current flow.
4. Apply the KneeStim and electrode.
5. Place the electrode pad around the same limb as the KneeStim. Positioning of this electrode will affect the sensation felt in the knee.
6. After use, return the electrode pad to the waxed paper sheet and seal it in the polythene bag supplied.

If necessary, hand wash the KneeStim and dry flat at room temperature.

CAUTION

If being used for pain management for Diabetic Neuropathy or other conditions where there is reduced or no sensory sensation then there is a small risk that reduced or uneven electrical contact may cause irritated or broken skin. Care needs to be taken to monitor the skin condition under the garment at frequent intervals.